b.safe

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www.bsafe.de
+49 (0) 93 46-92 86-0

Individual built-in systems
for safe solvent disposal
and waste management in
the laboratory

www.bsafe.de/en
+49 (0) 93 46-92 86-0

b.safe Displays

Version	No. of containers to be monitored	Dimensions	Cat-No.
Wall	6	154x53x99 mm	R 650-01
Wall	2	154x53x99 mm	R 652-01
Table	6	195x90x134 mm	R 655-01
Table	2	195x90x134 mm	R 657-01

b.safe Signal Cabels

Cable length	Connection 1	Connection 2	Cat-No.
2,0 m	socket	plug	R 645-02
3,0 m	socket	plug	R 645-03
5,0 m	socket	plug	R 645-05
10,0 m	socket	plug	R 645-10

Consumables & Accessories

b.safe Exhaust Filter

Size	Service life	Packing Unit	Cat.No.
Small	3 months	2 pcs.	M 506-01
Medium	6 months	2 pcs.	M 506-02
Large	12 months	2 pcs.	M 506-03

b.safe Mini-Display

Version	No. of containers to be monitored	Dimensions	Cat.No.
DLC (Wall)	1	84x84x70 mm	R 658-03
DLC (Table)	1	110x65x30 mm	R 659-03

6. As soon as all parts are dry, the glass/
the level indicator can be screwed
into the b-safe Waste Cap again.

4. Clean the protective glass / the level indicator in suitable baths (e.g. ultra sonic bath).

5. The protective glass of the b.safe Waste Cap LC can only be cleaned from below! Please use a small brush to get into the glass from below through the cap. **Attention:** don't use a metal brush!

3. Remove the protective glass / the level indicator from the b.safe Waste Cap LC/DLC
 - Attention:**
 - Don't immerse the b.safe Waste Cap with signal cable into liquids!
 - Clean the b.safe Waste Caps with a wet cloth
 - Do not disassemble the protective glass any further
 - Do not remove the signal cable from the b.safe Waste Cap


Cleaning & Maintenance

1. Possible labelling or adhesive residues on b.safe Waste Caps can be removed with isopropyl.
2. Once the b.safe Service Life Indicator is activated, the white inspection window turns into red with increasing time of use. If the white writing "CHANGE" is completely visible on the red background, the b.safe Exhaust filter is dead and has to be replaced.
1. Clean the b.safe Waste Caps with a wet cloth and, if applicable, with a rinsing agent and warm water.
Attention: the b.safe Exhaust Filter must not get wet. We recommend to remove it before cleaning
2. Unscrew the floater from the Waste Caps LC/DLC regularly and clean it by putting it into warm water (if necessary with rinsing agent). It is enough to flush it with clear water afterwards before mounting it again.

Connection of b.safe Waste Caps DLC

Connection of b.safe Waste Caps LC

Connection of b.safe Waste Caps with Pre- and Main Alarm


2. Connect signal cables of the b.safe Isolated Amplifier with the b.safe Display. Make sure not to mix up pre- and main alarm cables.

1. Connect signal cables of the b.safe Waste Cap with the b.safe Waste Cap with Isolated Amplifier.

EX

3. Connect the power supply with the power frid.

Upper row: main alarm (MAIN)
Lower row: pre alarm (PRE)


Important: Use b.safe Blind Fittings for gas tight closing of unused openings on a b.safe Waste Cap!

Before using Waste Caps LC/DLC always check if the application of an Isolated Amplifier (Cat. No. R 660-01) is necessary! Only use b.safe Signal Cables for the connection of Waste Caps LC/DLC and make sure that the plugs are not damaged (first plug and then screw)

Different kinds of connection of Waste Caps LC/DLC are explained in the instruction manuals of our Displays (e.g. Cat. No. R 657-01) or Mini-Displays (e.g. R 659-03).


b.safe Isolated Amplifier for galvanic isolation of the sensor side and the output side. Ignitable sparks caused by the components of the alerting systems are prevented in the explosive atmosphere on sensor side.


8. If you assemble b:safe System Waste Caps DLC 560 (Cat. No. R 1356-20) or b:safe System Waste Caps LC 560 (Cat. No. R 1306-01) to the b:safe Canisters 560/61 with inclined necks (Cat. No. R 236-10 and R 236-20) please consider that the level indicator always points upwards!

4. Tighten the b.safe Waste Cap onto a suitable canister until the insert does not move freely any longer.


5. Screw a b.safe Exhaust Filter into the free thread connection of the b.safe Waste Cap.

6. Push the button to activate the Service Life Indicator.
Optional: Note the exchange date on the back side.


7. Mark the b.safe Waste Cap with a foil pen or an adhesive label if necessary.

1. Push the capillaries through a suitable b.safe Fitting and afterwards through the b.safe Waste Cap.

2. Push the capillaries through the b.safe Fittings until their ends protrude on the bottom side of the b.safe Waste Cap for about 10 mm. Tighten the b.safe Fittings.

3. Screw the b.safe Hose Connector into the suitable thread of the b.safe Waste Cap. Push the tubing on the hose connector until it sticks.

Content

Functional instructions	3
Connection	3
Cleaning and Maintenance	6
Cleaning	6
More b.safe SYSTEM products	7

b.safe

Installation Manual b.safe SYSTEM Waste Caps LC/DLC