

## Model: ACID LINE

### Order number: 20-200667-031

Correct use: Storage non-flammable media



Classical door technology (door handle) door hung right For the storage of non-flammable liquids  
Type testing by MPA/TÜV SÜD according to DIN EN 16121 and DIN EN 16122, AfPS GS 2019: 01  
PAK Mark of approval: GS/CE mark, high quality mark of approval by TÜV SÜD Construction and  
colour:

- carcass of decor panels, RAL 7035
- door of decor panel, RAL 7035
- Fittings and features:
- tested object hinges (DIN EN 15570)
- cylinder lock for safe access
- 6 pull-out trays, PE
- exhaust air connection on the cabinet rear
- adjustable feet for level adjustment
- set back base, RAL 7015
- Further technical information:
- [www.dueperthal.com](http://www.dueperthal.com)

|                   |                          |
|-------------------|--------------------------|
| Category:         | Non-flammable substances |
| Model size        | M                        |
| Model type        | C PRO                    |
| Outer width [mm]: | 600                      |
| Outer depth [mm]: | 609                      |

|                                           |                      |
|-------------------------------------------|----------------------|
| Outer height [mm]:                        | 2.045                |
| Net weight [kg]:                          | 115,00               |
| Door finish                               | one piece            |
| Number of storage levels                  | 6                    |
| Number of storage areas per storage level | 1                    |
| Door operation                            | classical technology |
| Door technology                           | wing door technology |