

SPECIFICATIONS

Handle: Metal, BROEN-LAB POLYCOAT powder coated with media indication according to EN13792

Hose nozzle: Removable plastic

Headwork: Diaphragm type

Opening/closing function: 1 1/2 x 360°

Temperature range: 0-65°C

Max. working pressure: 10 bar / 147 psi

Test pressure: 1 x working pressure

Weight: 0.7 kg

MATERIALS WITH MEDIA CONTACT

Diaphragm: PTFE

Plastics: PVDF

Rubber: EPDM



ORDERING INFORMATION

For standard configuration please refer to the item numbers below.



Other options:

- different sizes of push fittings
- different inlet lengths

Model: With R3/8 push fitting

Media		S = 150mm	S = 200mm
 Distilled water		2533 046 1003-09	2533 033 1003-09
 Deionized water, cold		2533 046 1036-09	2533 033 1036-09

Model: Without push fitting

Media		S = 150mm	S = 200mm
 Distilled water		2533 008 1003-09	2533 010 1003-09
 Deionized water, cold		2533 008 1036-09	2533 010 1036-09