


ONE HANDLE MIXER WITH SWIVEL SPOUT



SPECIFICATIONS

Handle: Metal, BROEN-LAB POLYCOAT powder coated

Spout: 110° limited swivel action

Hose nozzle: Removable polypropylene (aerator included)

Cartridge: Ceramic (temperature and flow regulation type)

Opening/closing function: 25°, mixing function: ±45°

Temperature range: 0°C - 65°C

Max. working pressure: 10 bar / 147 psi

Weight: 2.8 kg

MATERIALS WITH MEDIA CONTACT

Metals: Brass

Plastics: PP, POM, PTFE, SPX (hose)

Sealing: EPDM

Lubricant: Silicone based

Ceramics: Al₂O₃ 96%

ORDERING INFORMATION

For standard configuration please refer to the item numbers below.

Other options:

- no swivel action
- full swivel action

Media:



Model

S = 200mm



Long Handle

G350318135

For further information and mounting instructions please refer to www.BROEN-LAB.com or contact your local BROEN-LAB representative.