



# ONE HANDLE TWO HOLE MIXER WITH SWIVEL SPOUT



## SPECIFICATIONS

**Handle:** Metal, BROEN-LAB POLYCOAT powder coated

**Spout:** Swivel action

**Hose nozzle:** Aerator

**Cartridge:** Ceramic (temperature and flow limiting type)

**Opening/closing function:** 25°, mixing function: ±45°

**Temperature range:** 0°C - 90°C

**Max. working pressure:** 10 bar / 147 psi

**Test pressure:** 1 x working pressure

**Inlets:** Eccentric couplings with male G1/2

**Weight:** 1.5 kg

## MATERIALS WITH MEDIA CONTACT

**Metals:** Brass

**Plastics:** POM, PTFE, SPX (hose)

**Sealing:** EPDM

**Lubricant:** Silicone based

**Ceramics:** Al<sub>2</sub>O<sub>3</sub> 96%

## ORDERING INFORMATION



For standard configuration please refer to the item numbers below.

Other options:

- removable hose nozzle
- bench mounted model
- different inlet connections  
(with BROEN-LAB UniFlex™ rear wall connection, lead-in, etc)

Media:



Model	S = 200mm	S = 250mm	S = 300mm
 Short Handle	2535 010 1035	2535 004 1035	2535 016 1035
 Long Handle	2535 007 1035	2535 001 1035	2535 013 1035

For further information and mounting instructions please refer to [www.BROEN-LAB.com](http://www.BROEN-LAB.com) or contact your local BROEN-LAB representative.