


SINGLE HANDLE MIXER BENCH MOUNTED



SPECIFICATIONS

- Handle:** Metal, BROEN-LAB POLYCOAT powder coated  
**Spout:** Swivel action  
**Hose nozzle:** Aerator  
**Cartridge:** Ceramic (temperature and flow regulation type)  
**Opening/closing function:** 25°, mixing function: ±45°  
**Temperature range:** 0°C - 90°C  
**Max. working pressure:** 10 bar / 147 psi  
**Test pressure:** 1 x working pressure  
**Weight:** 2.4 kg

MATERIALS WITH MEDIA CONTACT

- Metals:** Brass  
**Plastics:** POM, PTFE, SPX (hose)  
**Sealing:** EPDM  
**Lubricant:** Silicone based  
**Ceramics:** Al<sub>2</sub>O<sub>3</sub> 96%

ORDERING INFORMATION

For standard configuration please refer to the item numbers below.

Other options:

- integrated non-return valve
- removable and fixed hose nozzles
- flow restrictor (5 l/min)
- alternative hose inlet connections
- OD 10mm Cu pipes

Media:



Model	S = 150mm	S = 200mm	S = 250mm	S = 300mm
	2535 038 1035-52B	2535 039 1035-52B	2535 040 1035-52B	2535 041 1035-52B

For further information and mounting instructions please refer to [www.BROEN-LAB.com](http://www.BROEN-LAB.com) or contact your local BROEN-LAB representative.