


TURRET FITTINGS

With fixed nozzle



SPECIFICATIONS

Handle: Polypropylene with media indication according to EN13792

Hose nozzle: Metal, fixed (can be unscrewed with Allen key, thread: male G3/8)

Headwork: Needle

Opening/closing function: 3 x 360°

Max. working pressure:

16 bar / 232 psi

Allowable test pressure after installation:

1.5 x max working pressure without function of the valve

MATERIALS WITH MEDIA CONTACT

Materials with media contact are Oxygen approved and ultra-sonically cleaned.

Metals: Brass

Plastics: PVDF, PTFE (hose)

Sealing: FKM

Lubricant: Perfluoropolyether based


ORDERING INFORMATION

For standard configuration please refer to the item numbers below.

Other options:

- different media
- removable hose nozzle with union nut
- disability-friendly handle (single valve only)
- micro flow headwork
- alternative hose inlet connection / without a hose
- other inlet options (type, length)

With quick release coupling



TURRET FITTINGS

Model: with fixed nozzle

Media	 1-valve	 2-valve 90°	 2-valve 180°	 4-valve
 Compressed air	2532 075 1021-52B	2532 063 1021-52B	2532 076 1021-52B	2532 091 1021-52B
 Oxygen	2532 075 1022-52B	2532 063 1022-52B	2532 076 1022-52B	2532 091 1022-52B
 Nitrogen	2532 075 1023-52B	2532 063 1023-52B	2532 076 1023-52B	2532 091 1023-52B
 Carbon dioxide	2532 075 1024-52B	2532 063 1024-52B	2532 076 1024-52B	2532 091 1024-52B
 Argon	2532 075 1025-52B	2532 063 1025-52B	2532 076 1025-52B	2532 091 1025-52B
 Helium	2532 075 1026-52B	2532 063 1026-52B	2532 076 1026-52B	2532 091 1026-52B
 Dinitrogen monoxide	2532 075 1027-52B	2532 063 1027-52B	2532 076 1027-52B	2532 091 1027-52B

Fittings are delivered with a 700mm PTFE hose with OD10mm inlet connection.

For other options please refer to page 18 or contact your local BROEN-LAB representative.

Model: with quick release coupling (hexagon coded)

Media	 1-valve	 2-valve 90°	 2-valve 180°	 4-valve
 Compressed air	2532 127 1021-52B	2532 128 1021-52B	2532 129 1021-52B	2532 130 1021-52B
 Oxygen	2532 127 1022-52B	2532 128 1022-52B	2532 129 1022-52B	2532 130 1022-52B
 Nitrogen	2532 127 1023-52B	2532 128 1023-52B	2532 129 1023-52B	2532 130 1023-52B
 Carbon dioxide	2532 127 1024-52B	2532 128 1024-52B	2532 129 1024-52B	2532 130 1024-52B
 Argon	2532 127 1025-52B	2532 128 1025-52B	2532 129 1025-52B	2532 130 1025-52B
 Helium	2532 127 1026-52B	2532 128 1026-52B	2532 129 1026-52B	2532 130 1026-52B
 Dinitrogen monoxide	2532 127 1027-52B	2532 128 1027-52B	2532 129 1027-52B	2532 130 1027-52B

Fittings are delivered with a 700mm PTFE hose with OD10mm inlet connection.

For other options please refer to page 18 or contact your local BROEN-LAB representative.