


With fixed nozzle



With quick release coupling



SPECIFICATIONS

Handle: Polypropylene with media indication according to EN13792

Hose nozzle: Metal, fixed

Headwork: Ceramic, press/turn

Opening/closing function: 90°

Max. working pressure: 7 bar / 102 psi

Allowable test pressure after installation:
10 bar without function of the valve

MATERIALS WITH MEDIA CONTACT

Metals: Brass

Sealing: Nitrile (NBR)

Lubricant: Mineral oil based

Ceramics: Al₂O₃ 96%



ORDERING INFORMATION

For standard configuration please refer to the item numbers below.

Other options:

- different media
- flashback arrestor
- needle or 1/4-turn headwork
- disability-friendly handle (for needle and 1/4-turn headworks only)
- other inlet options (compression ring fittings, flexible SS pipe)
- national approvals

Models:

Media	with fixed nozzle	with quick release coupling
 Natural gas	2532 171 1031-52B	2532 203 1031-52B
 Liquefied petroleum gas	2532 171 1032-52B	2532 203 1032-52B

Fittings are delivered with a 700mm rubber hose with OD10mm inlet connection.

For other options please refer to page 19 or contact your local BROEN-LAB representative.