



## TWO VALVE BENCH MOUNTED FITTINGS ON A U-SHAPED COLUMN



700mm hose  
with OD10mm  
inlet connection

With wrist action levers  
(shown in CLOSED position)



**Model:** Polypropylene handle



### SPECIFICATIONS

**Handle:** Polypropylene, with media indication according to EN13792

**Hose nozzle:** Fixed Polypropylene (can be unscrewed by hand, thread: male G3/8)

**Headwork:** Rubber (wrist action lever: ceramic)

**Opening/closing function:** 2 x 360° (wrist action lever: 90°)

**Temperature range:** 0-65°C

**Max. working pressure:** 10 bar / 147 psi

**Test pressure:** 1 x working pressure

**Weight:** 2.5 kg

### MATERIALS WITH MEDIA CONTACT

**Metals:** Brass

**Plastics:** PP, SPX (hose)

**Sealing:** EPDM

**Lubricant:** Silicone based

**Ceramics:** Al<sub>2</sub>O<sub>3</sub> 96% (wrist action lever only)

### ORDERING INFORMATION

For standard configuration please refer to the item numbers below.

Other options:



- aerator
- removable polypropylene hose nozzle
- integrated non-return valve
- alternative hose inlet connection / without a hose
- other inlet options (type, length)
- national approvals
- WPH (Water Potable Hot)
- WNH (Water Non-Potable Hot)

BROEN-LAB UNIFLEX™ FITTINGS FOR WATER

Media	S = 150mm	S = 200mm
 Potable water, cold	2533 086 1001-52B	2533 149 1001-52B
 Non-potable water, cold	2533 086 1007-52B	2533 149 1007-52B

**Model:** Wrist action lever



Media	S = 150mm	S = 200mm
 Potable water, cold	2533 148 1001-52B	2533 150 1001-52B
 Non-potable water, cold	2533 148 1007-52B	2533 150 1007-52B

For further information and mounting instructions please refer to [www.BROEN-LAB.com](http://www.BROEN-LAB.com) or contact your local BROEN-LAB representative.