

## Model: UTS ERGO LINE

### Order number: 29-061467-005

Correct use: Passive storage



One-hand and classical door technology (door handle). For the storage of flammable liquids. Type testing by MPA/TÜV SÜD according to DIN EN 14470-1, DIN EN 16121 and DIN EN 16122, EK5/AK4 09-10:2009 and AfPS GS 2019: 01 PAK Mark of approval: GS/CE mark, high quality mark of approval by TÜV SÜD Construction and colour:

- Outer carcass sheet steel, RAL 7035
- Doors sheet steel, RAL 1018
- Inner carcass of decor panels, RAL 7035 Fittings and features:
- Doors close automatically in case of fire
- Doors remain open in any position
- Cylinder lock for safe access
- Centre partition wall
- 2x 1 automatically extending pull-out tray sheet steel, RAL 7035
- 2x 1 perforated sheet insert, RAL 7035
- Exhaust air connection socket at the rear of the cabinet
- Visual inspection of the ventilation opening
- Earthing connection at the rear of the cabinet
- Interior fittings conductively connected with carcass
- Classical slotted base, RAL 7015 Further technical information: [www.dueperthal.com](http://www.dueperthal.com)

|                    |                             |
|--------------------|-----------------------------|
| Category:          | Type 90 - flammable liquids |
| Model size         | XL                          |
| Outer width [mm]:  | 1.400                       |
| Outer depth [mm]:  | 594                         |
| Outer height [mm]: | 631                         |
| Net weight [kg]:   | 232,00                      |

|   |                      |
|---|----------------------|
| Number of storage levels                  | 1                    |
| Number of storage areas per storage level | 2                    |
| Door operation                            | classical technology |
| Door technology                           | wing door technology |